

Isothermal Master Mix - 300rxns

For research purposes only. Not for use in diagnostic procedures for clinical purposes.
FOR *IN VITRO* USE ONLY

Cat. No. ISO-DR004

Product Description:

The Isothermal Master Mix is an optimised Master Mix designed to simplify the preparation of an isothermal amplification assay. The master mix contains a proprietary fast, novel GspSSD2.0 DNA polymerase, proprietary thermostable inorganic pyrophosphatase, optimised reaction buffer, MgSO₄, dNTPs and a ds-DNA binding dye (FAM detection channel).

Stability:

The Isothermal Master Mix dried reagents are stable for a minimum of twelve months when stored at room temperature. Ideally store cold.

Supplied Reagents:

2 bottles, each containing 150rxn dried mastermix
1 bottle, each containing 5ml of resuspension buffer

Usage:

Add 2.25ml of the resuspension buffer to a dried mastermix bottle. Vortex gently to dissolve and leave to stand for 10mins before use.

A typical 25µl reaction would contain the following:

Component	Volume	Final concentration
Isothermal Master Mix	15.0µl	1x
Primers	5.0µl	-
Template	5.0µl	-
Total	25.0µl	

OptiGene Ltd.

Typical amplification conditions for this kit are:

Amplification	20mins	65°C
Anneal	98°C - 80°C - ramping at 0.05°C per second	

Note:

Due to the presence of the inorganic pyrophosphatase there will be no accumulation of PPI.

We recommend the following Primer levels per 25µl reaction

F3 & B3	5pmol each primer
LoopF & LoopB	10pmol each primer
FIP & BIP	20pmol each primer

Quality Control:

The Isothermal Master Mix has been tested for performance in Loop-mediated Isothermal DNA amplification.

Manufactured in the UK

OptiGene Ltd.
Unit 5 Blatchford Road
Horsham
WestSussex
RH13 5QR
UK

Tel: +44 (0) 1403 274980

FAX: +44 (0) 1403 271017

Notice to purchaser: limited license

The purchase price of this product includes a limited, non-transferable license under European Patents EP1,020,534, EP1,873,260, EP2,045,337, EP1,416,055, EP1,724,361 or their foreign counterparts, owned by Eiken Chemical Co., Ltd.