

# GspM2.0 DNA Polymerase

For research purposes only. Not for use in diagnostic procedures for clinical purposes. FOR *IN VITRO* USE ONLY

<b>Cat. No.</b>	<b>GSPM2.0-002</b>	<b>8000u</b>
<b>Cat. No.</b>	<b>GSPM2.0-003</b>	<b>40000u</b>

---

## Product Description

A recombinant engineered large fragment engineered DNA polymerase from *Geobacillus species M* expressed in *E.coli*.

**Applications:** This enzyme is designed for use in isothermal DNA amplification.

**Stability:** The undiluted enzyme mix is stable when stored at -20°C in a constant temperature freezer.

## Supplied Reagents:

Enzyme @ 8u/μl

10x Buffer (200mM Tris-HCl (pH 8.8 @ 25°C), 100mM KCl, 100mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> and 0.1% Triton X-100

5M Betaine

50mM MgSO<sub>4</sub>·7H<sub>2</sub>O.

**Storage Buffer:** 10mM Tris-HCl (pH 7.5), 50mM KCl, 0.1mM EDTA, 1mM DTT, 0.1% Triton-X-100 (v/v) and 50% glycerol.

**Quality Control:** Each lot is tested for activity in an Isothermal DNA amplification reaction.

For maximal activity and specificity in an isothermal amplification reaction we advise using Betaine at a final concentration of 1M.

---

**OptiGene Ltd.**